

Product Information

Product SKU:	CAB17946	Gene ID:	84315	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	73kDa	Conjugate:	-
Calculated MW:	73kDa	Isotype:	IgG

Immunogen Information

Background: Enables guanyl-nucleotide exchange factor activity. Predicted to be involved in protein secretion. Predicted to act upstream of or within cellular iron ion homeostasis and protein transport. Part of Mon1-Ccz1 complex.

Recommended Dilution:	WB, 1:500 - 1:2000
Synonyms:	SAND1; MON1A
Purification Method:	Affinity purification
Immunogen:	Recombinant protein of human MON1A.
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.