

PSAD Rabbit Polyclonal Antibody



CAB17958

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Protein Background

Immunogen information

Gene ID:

953220

Uniprot

Applications:

WB

Synonyms:

PSAD

Reactivity:

Synechocystis sp. PCC 6803

Immunogen:

A synthetic peptide of Synechocystis sp. PCC 6803 PSAD.

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Product Images
