

CAB17961

Product Information

Product SKU:	CAB17961	Gene ID:	73523	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	26kDa	Conjugate:	-
Calculated MW:	27kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be located in extracellular region. Is expressed in bone and rib. Orthologous to human PEBP4 (phosphatidylethanolamine binding protein 4).

Recommended Dilution:	WB,1:500 - 1:2000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200
Synonyms:	PEBP-4; PEBP4
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 90-242 of mouse PEBP4 (NP_082836.2).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.