

## CAB17963

---

**Product Information**

<b>Product SKU:</b>	CAB17963	<b>Gene ID:</b>	5900	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

**Additional Information**

<b>Observed MW:</b>	110kDa	<b>Conjugate:</b>	-
<b>Calculated MW:</b>	101kDa	<b>Isotype:</b>	IgG

---

**Immunogen Information**

**Background:** Guanine nucleotide dissociation stimulators (GDSs, or exchange factors), such as RALGDS, are effectors of Ras-related GTPases (see MIM 190020) that participate in signaling for a variety of cellular processes.

**Recommended Dilution:** WB,1:500 - 1:2000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

**Synonyms:** RGF; RGDS; RalGEF; RALGDS

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 771-902 of human RALGDS (NP\_001258705).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.