

CAB17966

Product Information

Product SKU:	CAB17966	Gene ID:	4519	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	43kDa	Conjugate:	-
Calculated MW:	43kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable metal ion binding activity. Predicted to be involved in several processes, including electron transport coupled proton transport; response to cobalamin; and response to glucagon. Located in mitochondrion. Implicated in ovarian carcinoma and urinary bladder cancer.

Recommended Dilution: WB,1:500 - 1:2000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:100

Synonyms: MTCYB; CYTB; MT-CYB

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human MT-CYB (YP_003024038.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.