

MT-ND5 Rabbit Polyclonal Antibody



CAB17972

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

67kDa

Calculated MW:

67kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

4540

Uniprot

P03915

Synonyms:

MTND5; MT-ND5


Immunogen:

A synthetic peptide corresponding to a sequence within amino acids 510 to the C-terminus of human MT-ND5 (BAG14520.1).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - MT-ND5 Rabbit pAb (CAB17972)