

CAB17978

Product Information

Product SKU:	CAB17978	Gene ID:	2308/2309/4303	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	80kDa	Conjugate:	-
Calculated MW:	-	Isotype:	IgG

Immunogen Information

Background: The Forkhead family of transcription factors mediates many aspects of physiology, including stress response, metabolism, commitment to apoptosis, and development. The Forkhead box subfamily O (FoxO) proteins have garnered particular interest due to their involvement in the modulation of cancer, metabolism and cardiovascular biology. In mammals, four subfamily members have been identified: FOXO1 (previously also known as FKHR), FOXO3 (also known as FKHL1), FOXO4 (also known as AFX), and FOXO6.

Recommended Dilution: WB, 1:500 - 1:2000 IHC-P, 1:50 - 1:200

Synonyms: -

Purification Method: Affinity purification

Immunogen: A synthetic peptide of human FOXO1/FOXO3A/FOXO4 (NP_002006.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.