

ZDHHC7 Rabbit Polyclonal Antibody



CAB17981

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

33kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

55625

Uniprot

Q9NXF8

Synonyms:

DHHC7; SERZ-B; SERZ1; ZNF370; ZDHHC7


Immunogen:

Recombinant protein of human ZDHHC7.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - ZDHHC7 Rabbit pAb (CAB17981)