

# ZDHHC20 Rabbit Polyclonal Antibody



CAB17982

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

42kDa

### Calculated MW:

### Applications:

WB

### Reactivity:

Human, Mouse, Rat

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

253832

#### Uniprot

Q5W0Z9

#### Synonyms:

4933421L13Rik; DHHC-20; DHHC20; ZDHHC20

#### Immunogen:


Recombinant protein of human ZDHHC20.

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot - ZDHHC20 Rabbit pAb (CAB17982)