

MonoMethyl-eEF1A-K55 Rabbit Polyclonal Antibody

CAB17986



Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

Applications:

Reactivity:

Antibody Information

Recommended dilutions:

Source:

Rabbit

Isotype:

N/A

Purification:

N/A

Protein Background

Immunogen information

Gene ID:

Uniprot

Synonyms:

Immunogen:

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
