

Product Information

Product SKU:	CAB17987	Gene ID:	-	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	50kDa	Conjugate:	-
Calculated MW:	-	Isotype:	IgG

Immunogen Information

Background:

Recommended Dilution:	WB, 1:500 - 1:1000
Synonyms:	-
Purification Method:	Affinity purification
Immunogen:	A synthetic dimethylated peptide around K55 of human eEF1A (NP_001393.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.