

C10orf67 Rabbit Polyclonal Antibody



CAB17993

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

64kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

256815

Uniprot

Q8IYJ2

Synonyms:

C10orf115; LINC01552; bA215C7.4; C10orf67


Immunogen:

Recombinant protein of human C10orf67.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - C10orf67 Rabbit pAb (CAB17993)