

SENP2 Rabbit Polyclonal Antibody



CAB17994

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

68kDa

Calculated MW:

68kDa

Applications:

WB

Reactivity:

Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

59343

Uniprot

Q9HC62

Synonyms:

AXAM2; SMT3IP2; SENP2


Immunogen:

Recombinant protein of human SENP2.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - SENP2 Polyclonal Antibody (CAB17994)