

LAPTM5 Rabbit Polyclonal Antibody



CAB17995

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

30kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

7805

Uniprot

Q13571

Synonyms:

CLAST6; LAPTM5


Immunogen:

Recombinant protein of human LAPTM5.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - LAPT5 Polyclonal Antibody (CAB17995)