

Product Information

Product SKU:	CAB17995	Gene ID:	7805	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	25kDa	Conjugate:	-
Calculated MW:	30kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a transmembrane receptor that is associated with lysosomes. The encoded protein, also known as E3 protein, may play a role in hematopoiesis.

Recommended Dilution:	WB,1:500 - 1:1000
Synonyms:	CLAST6; LPTM5
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 153-262 of human LPTM5 (NP_006753.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.