

CAB18026

Product Information

Product SKU:	CAB18026	Gene ID:	8667	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	40kDa	Conjugate:	Unconjugated
Calculated MW:	40kDa	Isotype:	IgG

Immunogen Information

Background: Enables deubiquitinase activity. Contributes to translation initiation factor activity. Involved in negative regulation of proteasomal ubiquitin-dependent protein catabolic process and translational initiation. Located in extracellular exosome and membrane. Part of eukaryotic translation initiation factor 3 complex. Implicated in breast cancer; prostate cancer; and prostate carcinoma. Biomarker of prostate cancer.

Recommended Dilution: WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200

Synonyms: EIF3S3; eIF3-p40; eIF3-gamma; 3H

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 80-340 of human EIF3H (NP_003747.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.