

**[KO Validated] HERC4 Rabbit Polyclonal
Antibody**
CAB18049



Product Information
Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

118kDa

Calculated MW:

12kDa/13kDa/106kDa/109kDa/117kDa/118kDa

Applications:

WB

Reactivity:

Human, Mouse, Rat

Protein Background

Immunogen information

Gene ID:

26091

Uniprot

Q5GLZ8

Synonyms:

HERC4

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 400-700 of human HERC4 (NP_056416.2).

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG


Purification:

Affinity purification

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, pH7.3.

Product Images



Western blot analysis of extracts from normal (control) and HERC4 knockout (KO) HeLa cells, using HERC4 antibody (CAB18049) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (CABS014) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (CABM00020). Exposure time: 1min.