

CAB18060

Product Information

| | | | | | |
|---------------------|----------|----------------------|--------|--------------------|-----------------|
| Product SKU: | CAB18060 | Gene ID: | 7532 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Human,Mouse,Rat |

Additional Information

| | | | |
|-----------------------|-------|-------------------|--------------|
| Observed MW: | 28kDa | Conjugate: | Unconjugated |
| Calculated MW: | 28kDa | Isotype: | IgG |

Immunogen Information

Background: This gene product belongs to the 14-3-3 family of proteins which mediate signal transduction by binding to phosphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals, and this protein is 100% identical to the rat ortholog. It is induced by growth factors in human vascular smooth muscle cells, and is also highly expressed in skeletal and heart muscles, suggesting an important role for this protein in muscle tissue. It has been shown to interact with RAF1 and protein kinase C, proteins involved in various signal transduction pathways.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: DEE56; EIEE56; PPP1R170; 14-3-3GAMMA; ma

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human 14-3-3 gamma (P61981).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.