

CAB18101

Product Information

Product SKU:	CAB18101	Gene ID:	11168	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	75kDa	Conjugate:	Unconjugated
Calculated MW:	60kDa	Isotype:	IgG

Immunogen Information

Background: Enables DNA-binding transcription factor binding activity; chromatin binding activity; and transcription coactivator activity. Involved in mRNA 5'-splice site recognition and positive regulation of transcription by RNA polymerase II. Located in heterochromatin; nuclear periphery; and nucleoplasm.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: p52; p75; PAIP; DFS70; LEDGF; PSIP2; P1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 301-530 of human PSIP1 (NP_150091.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.