

LUC7L3 Rabbit Polyclonal Antibody



CAB18104

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

51747

Uniprot

O95232

Synonyms:

CRA; CREAP-1; CROP; LUC7A; OA48-18; hLuc7A; LUC7L3

Immunogen:

Recombinant fusion protein containing a sequence corresponding to amino acids 60-200 of human LUC7L3 (NP_057508.2).

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
