

K6-linkage Specific Polyubiquitin Rabbit Polyclonal Antibody

CAB18106



Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:**Calculated MW:****Applications:**

WB

Reactivity:

ALL

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

Uniprot


Synonyms:**Immunogen:**

A synthetic peptide of K6-linkageSpecificPolyubiquitin.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - K6-linkage Specific Polyubiquitin Rabbit pAb (CAB18106)