

---

## Product Information

<b>Product SKU:</b>	CAB18108	<b>Gene ID:</b>	10683	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

## Additional Information

<b>Observed MW:</b>	65kDa	<b>Conjugate:</b>	-
<b>Calculated MW:</b>	65kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** This gene encodes a member of the delta protein ligand family. This family functions as Notch ligands that are characterized by a DSL domain, EGF repeats, and a transmembrane domain. Mutations in this gene cause autosomal recessive spondylocostal dysostosis 1. Two transcript variants encoding distinct isoforms have been identified for this gene.

**Recommended Dilution:** WB,1:100 - 1:500

**Synonyms:** SCDO1; DLL3

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 520-618 of human DLL3 (NP\_058637.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.