
Product Information

Product SKU:	CAB18111	Gene ID:	2847	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	46kDa	Conjugate:	-
Calculated MW:	46kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene, a member of the G protein-coupled receptor family 1, is an integral plasma membrane protein which binds melanin-concentrating hormone. The encoded protein can inhibit cAMP accumulation and stimulate intracellular calcium flux, and is probably involved in the neuronal regulation of food consumption. Although structurally similar to somatostatin receptors, this protein does not seem to bind somatostatin.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: SLC1; GPR24; MCH1R; SLC-1; MCH-1R; MCHR1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human MCHR1 (NP_005288.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.