

MCHR1 Rabbit Polyclonal Antibody



CAB18111

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

46kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

2847

Uniprot

Q99705

Synonyms:

GPR24; MCH-1R; MCH1R; SLC-1; SLC1; MCHR1


Immunogen:

Recombinant protein of human MCHR1.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - MCHR1 Polyclonal Antibody (CAB18111)