

ZDHHC6 Rabbit Polyclonal Antibody



CAB18113

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

56kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

64429

Uniprot

Q9H6R6

Synonyms:

DHHC-6; ZNF376; ZDHHC6


Immunogen:

Recombinant protein of human ZDHHC6.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - ZDHHC6 Rabbit pAb (CAB18113)