

---

## Product Information

<b>Product SKU:</b>	CAB18117	<b>Gene ID:</b>	5604/5605	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Mouse,Rat

---

## Additional Information

<b>Observed MW:</b>	44kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	40kDa/43kDa/44kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

<b>Background:</b>	The protein MEK1 and MEK2 are two members of the dual specificity protein kinase family, which act as mitogen-activated protein (MAP) kinase kinase. MAP kinases, also known as extracellular signal-regulated kinases (ERKs), act as an integration point for multiple biochemical signals. This protein kinase lies upstream of MAP kinases and stimulates the enzymatic activity of MAP kinases upon wide variety of extra- and intracellular signals. As an essential component of MAP kinase signal transduction pathway, this kinase is involved in many cellular processes such as proliferation, differentiation, transcription regulation and development.
<b>Recommended Dilution:</b>	WB,1:100 - 1:500
<b>Synonyms:</b>	MEK1/MEK2
<b>Purification Method:</b>	Affinity purification
<b>Immunogen:</b>	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human MAP2K1 (NP_002746.1).
<b>Storage:</b>	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.