

RXRB Rabbit Polyclonal Antibody



CAB18119

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

70kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

6257

Uniprot

P28702

Synonyms:

DAUDI6; H-2RIIBP; NR2B2; RCoR-1; RXRB


Immunogen:

Recombinant protein of human RXRB.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - RXRB Polyclonal Antibody (CAB18119)