

YTHDF1 Rabbit Polyclonal Antibody



CAB18126

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

61kDa

Calculated MW:

20kDa/60kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

54915

Uniprot

Q9BYJ9

Synonyms:

YTHDF1; C20orf21


Immunogen:

Recombinant protein of human YTHDF1.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - YTHDF1 Polyclonal Antibody (CAB18126)