

CAB18127

Product Information

Product SKU:	CAB18127	Gene ID:	2534	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	61kDa	Conjugate:	Unconjugated
Calculated MW:	61kDa	Isotype:	IgG

Immunogen Information

**Background:** This gene is a member of the protein-tyrosine kinase oncogene family. It encodes a membrane-associated tyrosine kinase that has been implicated in the control of cell growth. The protein associates with the p85 subunit of phosphatidylinositol 3-kinase and interacts with the fyn-binding protein. Alternatively spliced transcript variants encoding distinct isoforms exist.

**Recommended Dilution:** WB,1:100 - 1:500 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

**Synonyms:** SLK; SYN; p59-FYN; Fyn

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 200-300 of human Fyn (NP\_002028.1).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.