
Product Information

Product SKU:	CAB18133	Gene ID:	55908	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	24kDa	Conjugate:	-
Calculated MW:	22kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable hormone activity. Involved in regulation of lipid metabolic process and triglyceride homeostasis. Acts upstream of or within positive regulation of protein processing and regulation of lipoprotein metabolic process. Located in extracellular region.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: RIFL; TD26; PRO1185; PVPA599; C19orf80; ANGPTL8

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 22-198 of human ANGPTL8 (NP_061157.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.