

# TPST2 Rabbit Polyclonal Antibody



CAB18141

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

41kDa

### Applications:

WB

### Reactivity:

Mouse

## Protein Background

The protein encoded by this gene catalyzes the O-sulfation of tyrosine residues within acidic regions of proteins. The encoded protein is a type II integral membrane protein found in the Golgi body. Two transcript variants encoding the same protein have been found for this gene.

## Immunogen information

### Gene ID:

8459

### Uniprot

O60704

### Synonyms:

TPST2; TANGO13B

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:1000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

### Immunogen:

Recombinant protein of human TPST2

### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, pH7.3.

## Product Images

---