

RAB11FIP5 Rabbit Polyclonal Antibody



CAB18142

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

65kDa

Calculated MW:

70kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

26056

Uniprot

Q9BXF6

Synonyms:

RAB11FIP5; GAF1; RIP11; pp75


Immunogen:

Recombinant protein of human RAB11FIP5

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, pH7.3.

Product Images



Western blot - RAB11FIP5 Polyclonal Antibody (CAB18142)