

Product Information

Product SKU:	CAB18142	Gene ID:	26056	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	70-90kDa	Conjugate:	Unconjugated
Calculated MW:	70kDa	Isotype:	IgG

Immunogen Information

Background: Enables gamma-tubulin binding activity. Involved in cellular response to acidic pH; negative regulation of adiponectin secretion; and regulation of protein localization to cell surface. Located in centriolar satellite and mitochondrial outer membrane.

Recommended Dilution:	WB,1:500 - 1:1000 IF/ICC,1:50 - 1:200
Synonyms:	GAF1; pp75; RIP11; gaf-1; rab11-FIP5; RAB11FIP5
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 354-653 of human RAB11FIP5 (NP_056285.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.