

MUC5B Rabbit Polyclonal Antibody



CAB18146

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

Applications:

IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

727897

Uniprot

Q9HC84

Synonyms:

MG1; MUC-5B; MUC5; MUC9; MUC5B

Immunogen:

Recombinant protein of human MUC5B.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
