

AJUBA Rabbit Polyclonal Antibody



CAB18153

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

57kDa

Calculated MW:

13kDa/56kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

AJUBA (ajuba LIM protein) is a protein-coding gene. GO annotations related to this gene include chromatin binding and actin filament binding. An important paralog of this gene is LPP.

Immunogen information

Gene ID:

84962

Uniprot

Q96IF1

Synonyms:

AJUBA; JUB


Immunogen:

Recombinant Protein of human AJUBA.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - AJUBA Polyclonal Antibody (CAB18153)