

MYMK Rabbit Polyclonal Antibody



CAB18158

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

27KDa

Calculated MW:

25kDa

Applications:

WB

Reactivity:

Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

389827

Uniprot

A6NI61

Synonyms:

TMEM8C; TMEM226; MYOMAKER; MYMK

Immunogen:

A synthetic peptide of human MYMK.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot - MYMK Rabbit pAb (CAB18158)

