

OSR2 Rabbit Polyclonal Antibody



CAB18178

Product Information

Size:

50uL, 100uL, 200uL

Observed MW:

30kDa

Calculated MW:

Applications:

WB

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

116039

Uniprot

Q8N2R0

Synonyms:

OSR2

Immunogen:

Recombinant protein of human OSR2.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

Western blot - OSR2 Rabbit pAb (CAB18178)

