

## CAB18181

---

**Product Information**

<b>Product SKU:</b>	CAB18181	<b>Gene ID:</b>	5829	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse

---

**Additional Information**

<b>Observed MW:</b>	65kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	65kDa	<b>Isotype:</b>	IgG

---

**Immunogen Information**

**Background:** This gene encodes a cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extracellular matrix (focal adhesion). Alternatively spliced transcript variants encoding different isoforms have been described for this gene. These isoforms exhibit different expression pattern, and have different biochemical, as well as physiological properties (PMID:9054445).

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** Paxillin; PXN; paxillin

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Paxillin (NP\_002850.2).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.