

Mrpl46 Rabbit Polyclonal Antibody



CAB18187

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

32kDa

Calculated MW:

32kDa

Applications:

WB IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IHC 1:50
- 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

67308

Uniprot

Q9EQI8

Synonyms:

C80107; LIECG2; P2ECSL; C15orf4; 3110052F15Rik; Mrpl46


Immunogen:

Recombinant protein of human Mrpl46 .

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - Mrpl46 Rabbit pAb (CAB18187)