

K33-linkage Specific Polyubiquitin Rabbit Polyclonal Antibody

CAB18199



Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Calculated MW:

Applications:

WB

Reactivity:

ALL

Protein Background

Immunogen information

Gene ID:

Uniprot

Synonyms:

Immunogen:

Recombinant protein of human K33-linkage Specific Polyubiquitin.

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG


Purification:

Affinity purification

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - K33-linkage Specific Polyubiquitin Rabbit pAb (CAB18199)