

Mrpl58 Rabbit Polyclonal Antibody



CAB18209

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

23kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

68572

Uniprot

Q8R035

Synonyms:

Ict1; MRP-L58; 1110001A02Rik; 1110002E03Rik; Mrpl58

Immunogen:

Recombinant protein of human Mrpl58.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
