

IGF1R/INSR Rabbit Polyclonal Antibody



CAB18213

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

154kDa/155kDa/156kDa

Applications:

WB

Reactivity:

Human

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

3480/3643

Uniprot

P08069/P06213

Synonyms:

IGF1R/INSR

Immunogen:

Recombinant protein of human IGF1R/INSR.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
