

# EID1 Rabbit Polyclonal Antibody



**CAB18217**

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

21kDa

### Applications:

WB

### Reactivity:

Human

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

23741

#### Uniprot

Q9Y6B2

#### Synonyms:

CRI1; EID-1; RBP21; PTD014; C15orf3; PNAS-22; IRO45620; EID1

#### Immunogen:

Recombinant protein of human EID1.

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---