

CAB18218

Product Information

| | | | | | |
|--------------|----------|---------------|--------|-------------|-----------------|
| Product SKU: | CAB18218 | Gene ID: | 23035 | Size: | 20uL, 100uL |
| Clone No: | - | Host Species: | Rabbit | Reactivity: | Human,Mouse,Rat |

Additional Information

| | | | |
|----------------|--------|------------|-----|
| Observed MW: | 140kDa | Conjugate: | - |
| Calculated MW: | 147kDa | Isotype: | IgG |

Immunogen Information

Background: Predicted to enable protein serine/threonine phosphatase activity. Predicted to be involved in signal transduction. Located in several cellular components, including intercellular bridge; mitotic spindle; and nucleoplasm.

| | |
|-----------------------|---|
| Recommended Dilution: | WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200 |
| Synonyms: | PPM3B; PHLPPL; PHLPP2 |
| Purification Method: | Affinity purification |
| Immunogen: | A synthetic peptide corresponding to a sequence within amino acids 1250-1323 of human PHLPP2 (NP_055835.2). |
| Storage: | Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3. |