

PHLPP2 Rabbit Polyclonal Antibody



CAB18218

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

140kDa

Calculated MW:

147kDa

Applications:

WB IF

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000 IF 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

23035

Uniprot

Q6ZVD8

Synonyms:

PHLPP2; PHLPL; PPM3B


Immunogen:

Recombinant protein of human PHLPP2.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot - PHLPP2 Polyclonal Antibody (CAB18218)