

HILPDA Rabbit Polyclonal Antibody



CAB18221

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

7kDa

Applications:

IHC

Reactivity:

Human, Mouse, Rat

Antibody Information

Recommended dilutions:

IHC 1:50 - 1:200

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Protein Background

Immunogen information

Gene ID:

29923

Uniprot

Q9Y5L2

Synonyms:

HIG2; HIG-2; C7orf68; HILPDA


Immunogen:

Recombinant protein of human HILPDA.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry - HILPDA Polyclonal Antibody (CAB18221)