

Product Information

Product SKU:	CAB18221	Gene ID:	29923	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	Refer to figures	Conjugate:	-
Calculated MW:	7kDa	Isotype:	IgG

Immunogen Information

Background: Enables signaling receptor binding activity. Involved in several processes, including autocrine signaling; cellular response to hypoxia; and positive regulation of lipid storage. Located in several cellular components, including cell surface; lipid droplet; and secretory granule.

Recommended Dilution:	IHC-P,1:50 - 1:200
Synonyms:	HIG2; HIG-2; C7orf68; HILPDA
Purification Method:	Affinity purification
Immunogen:	Recombinant fusion protein containing a sequence corresponding to amino acids 1-63 of human HILPDA (NP_037464.1).
Storage:	Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.