

# SCGB3A1 Rabbit Polyclonal Antibody



CAB18228

---

## Product Information

### Size:

20uL, 50uL, 100uL, 200uL

### Observed MW:

Refer to figures

### Calculated MW:

10kDa

### Applications:

WB

### Reactivity:

Human

## Antibody Information

### Recommended dilutions:

WB 1:500 - 1:2000

### Source:

Rabbit

### Isotype:

IgG

### Purification:

Affinity purification

## Protein Background

### Immunogen information

#### Gene ID:

92304

#### Uniprot

Q96QR1

#### Synonyms:

HIN1; HIN-1; LU105; UGRP2; PhSP-2; SCGB3A1

#### Immunogen:

Recombinant protein of human SCGB3A1.

#### Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---