

CLTCL1 Rabbit Polyclonal Antibody



CAB18230

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

187kDa

Applications:

WB

Reactivity:

Human

Protein Background

This gene is a member of the clathrin heavy chain family and encodes a major protein of the polyhedral coat of coated pits and vesicles. Chromosomal aberrations involving this gene are associated with meningioma, DiGeorge syndrome, and velo-cardio-facial syndrome. Multiple transcript variants encoding different isoforms have been found for this gene.

Immunogen information

Gene ID:

8218

Uniprot

P53675

Synonyms:

CLTCL1; CHC22; CLH22; CLTCL; CLTD

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

Recombinant protein of human CLTCL1.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
