

RIOK1 Rabbit Polyclonal Antibody



CAB18234

Product Information

Size:

20uL, 50uL, 100uL, 200uL

Observed MW:

Refer to figures

Calculated MW:

65kDa

Applications:

WB

Reactivity:

Human

Protein Background

The protein encoded by this gene competes with pICln for inclusion in the protein arginine methyltransferase 5 complex. This complex targets substrates for dimethylation. The encoded protein is essential for the last steps in the maturation of 40S subunits.

Immunogen information

Gene ID:

83732

Uniprot

Q9BRS2

Synonyms:

RIO1; AD034; RRP10; bA288G3.1; RIOK1

Antibody Information

Recommended dilutions:

WB 1:500 - 1:2000

Source:

Rabbit

Isotype:

IgG

Purification:

Affinity purification

Immunogen:

Recombinant protein of human RIOK1.

Storage:

Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images
